

Rabbit Anti-Human Cd85I/J Monoclonal Antibody For Flow Cytometry - O75019/Q8Nhl6

Item Number: CM0012700

	Introduction -conjugated rabbit monoclonal antibody targeting human CD85i/j (LILRA1/LILRB1). Validated for flow cytometry, recognizing native human antigens. Ideal for immune cell profiling and receptor studies. Learn More
--	---

Parameter	Value
Product Name	Rabbit anti-Human CD85i/j mAb
Remarks/Alias	LIR6; CD85I; LIR-6; CD85J; ILT-2; ILT2; LIR-1; LIR1; MIR-7; MIR7; PIR-B; PIRB
Species	Human
Gene ID (Human)	11024/10859
Gene ID	11024/10859
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75019/Q8NHL6
Protein Weight	53kDa
Shipping	Ice bag