


Rabbit Anti-Human Cd85j/Ilt2 Monoclonal Antibody For Fc - Q8Nhl6

Item Number: CM0011150



Introduction

-conjugated rabbit monoclonal antibody targeting human CD85j/ILT2 (LILRB1) for flow cytometry. Recognizes a broad spectrum of MHC class I molecules; ideal for studying immune checkpoint regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD85j/ILT2 mAb
Remarks/Alias	ILT2; LIR1; MIR7; PIRB; CD85j; ILT-2; LIR-1; MIR-7; PIR-B
Species	Human
Gene ID (Human)	10859
Gene ID	10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-458 of human CD85j/ILT2 (NP_006660.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8NHL6
Protein Weight	49kDa/70kDa/71kDa
Shipping	Ice bag