


Rabbit Anti-Human Cd85J/Ilt2 Monoclonal Antibody For Fc - Q8Nhl6

Item Number: CM0009508



Introduction

-conjugated anti-CD85j/ILT2 rabbit monoclonal antibody for flow cytometry. Targets human LILRB1, an inhibitory MHC class I receptor expressed on immune cells, enabling immune regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD85j/ILT2 mAb
Remarks/Alias	ILT2, LIR1, MIR7, PIRB, CD85J, ILT-2, LIR-1, MIR-7, PIR-B
Species	Human
GeneID (Human)	10859
GeneID	10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-458 of human CD85j/ILT2 (NP_006660.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8NHL6
Protein Weight	49kDa/70kDa/71kDa
Shipping	Ice bag