


Rabbit Anti-Human Cd87/Plaur Monoclonal Antibody For Flow Cytometry - Q03405

Item Number: CM0005244



Introduction

-conjugated rabbit monoclonal antibody targeting human CD87/PLAUR (urokinase plasminogen activator receptor), validated for flow cytometry. Enables study of cell surface uPAR in proteolysis, adhesion, and migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD87/PLAUR mAb
Remarks/Alias	PLAUR; CD87; U-PAR; UPAR; URKR; plasminogen activator, urokinase receptor
Species	Human
GeneID (human)	5329
GeneID	5329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-303 of human CD87/PLAUR (NP_002650.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q03405
Protein Weight	31kDa/32kDa/36kDa
Shipping	Ice bag