


Rabbit Anti-Human Cd89 Monoclonal Antibody For Fc - P24071

Item Number: CM0012692



Introduction

-conjugated rabbit monoclonal antibody targeting human CD89 (FcalphaR).
Recombinant immunogen, validated for flow cytometry (FC). Cross-reacts with human. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD89 mAb
Remarks/Alias	CD89; FcalphaR; FcalphaRI; CTB-61M7.2
Species	Human
GeneID (Human)	2204
GeneID	2204
Immunogen	Recombinant protein of human CD89.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24071
Protein Weight	32kDa
Shipping	Ice bag