


Rabbit Anti-Human Cd89 Monoclonal Antibody For Fc - P24071

Item Number: CM0012544



Introduction

-conjugated rabbit monoclonal antibody targeting human CD89 (FcalphaRI) for flow cytometry. Validated for human samples. UniProt accession: P24071. Suitable for myeloid cell and immune receptor research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD89 mAb
Remarks/Alias	CD89; FcalphaR; FcalphaRI; CTB-61M7.2
Species	Human
Gene ID (Human)	2204
Gene ID	2204
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-227 of human CD89 (NP_001991.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24071
Protein Weight	32kDa
Shipping	Ice bag