


Rabbit Anti-Human Cd89 Monoclonal Antibody For Fc - P24071

Item Number: CM0004381



Introduction

-conjugated rabbit anti-human CD89 monoclonal antibody for flow cytometry detection of IgA Fc receptor (FCAR) on human immune cells. High specificity, recombinant immunogen.

[Learn More](#)

Specification	Details
Product Name	Rabbit anti-Human CD89 mAb
Remarks/Alias	CD89; FcαRI; FcαRI; CTB-61M7.2
Species	Human
Gene ID (Human)	2204
Gene ID	2204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-227 of human CD89 (NP_001991.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P24071
Protein Weight	32kDa
Shipping	Ice bag