

Rabbit Anti-Human Cd8B Monoclonal Antibody For Flow Cytometry - P10966

Item Number: CM0006081

Introduction

-conjugated rabbit monoclonal antibody targeting human CD8b, validated for flow cytometry. CD8b is a T-cell coreceptor involved in MHC class I recognition and cytotoxic T-cell activation. Ideal for immune profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD8b mAb
Remarks/Alias	LY3; P37; LEU2; LYT3; Ly-3; CD8B1; CD8beta
Species	Human
GeneID (Human)	926
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-170 of human CD8b (NP_004922.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10966
Protein Weight	24kDa
Shipping	Ice bag