


Rabbit Anti-Human Cd9 Monoclonal Antibody For Fc - P21926

Item Number: CM0006236



Introduction

Rabbit anti-Human CD9 monoclonal antibody validated for flow cytometry, targeting the 25 kDa tetraspanin CD9 (P21926). Ideal for immune cell and exosome studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD9 mAb
Remarks/Alias	MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29
Species	Human
Gene ID	928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human CD9 (NP_001760.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21926
Protein Weight	25kDa
Shipping	Ice bag