


Rabbit Anti-Human Cd9 Monoclonal Antibody () For Flow Cytometry - P21926

Item Number: CM0012670



Introduction

-conjugated rabbit monoclonal antibody against human CD9, validated for flow cytometry. Targets tetraspanin CD9 (TSPAN29) involved in cell adhesion, sperm-egg fusion, platelet activation, and tumor metastasis. For research use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD9 mAb
Remarks/Alias	MIC3, MRP-1, BTCC-1, DRAP-27, TSPAN29, TSPAN-29
Species	Human
Gene ID (Human)	928
Gene ID	928
Immunogen	Recombinant protein / Recombinant fusion protein containing amino acids 1-228 of human CD9 (NP_001760.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21926
Protein Weight	25kDa
Shipping	Ice bag