

Rabbit Anti-Human Cd90/Thy1 Mab For Fc - P04216

Item Number: CM0004204

Introduction

-conjugated rabbit monoclonal antibody targeting human CD90/Thy1 (Thy-1 membrane glycoprotein). Validated for flow cytometry. Useful for neuroscience and immunology research requiring detection of this GPI-anchored cell surface antigen.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD90/Thy1 mAb
Remarks/Alias	CD90; CDw90
Species	Human
GeneID (Human)	7070
GeneID	7070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CD90/Thy1 (NP_006279.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag