

Rabbit Anti-Human Cd90/Thy1 Monoclonal Antibody For Flow Cytometry - P04216

Item Number: CM0009636

Introduction

-conjugated rabbit monoclonal antibody against human CD90/Thy1 (Swiss-Prot P04216), validated for flow cytometry. Suitable for studies of cell membrane glycoprotein involved in synaptogenesis and cell interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD90/Thy1 mAb
Remarks/Alias	CD90, CDw90
Species	Human
GeneID (Human)	7070
GeneID	7070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CD90/Thy1 (NP_006279.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag