


Rabbit Anti-Human Cd93 Mab For Fc - Q9Npy3

Item Number: CM0009018



Introduction

Rabbit anti-Human CD93 monoclonal antibody for flow cytometry (FC) applications. Validated for human samples, this antibody targets the complement component C1q receptor, ideal for studying phagocytosis and endothelial cell function.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD93 mAb
Remarks/Alias	C1QR1; C1qRP; CDw93; ECSM3; MXRA4; C1qR(P); dj737E23.1
Species	Human
GeneID (Human)	22918
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-580 of human CD93 (NP_036204.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NPY3
Protein Weight	69kDa
Shipping	Ice bag