


Rabbit Anti-Human Cd93 Monoclonal Antibody For Fc - Q9Npy3

Item Number: CM0009729



Introduction

This rabbit monoclonal antibody against human CD93, conjugated to , is validated for flow cytometry. Ideal for studying CD93 expression in endothelial and myeloid cells. Predicted molecular weight 69 kDa (UniProt Q9NPY3).

[Learn More](#)

Parameter	Description
Product Name	Rabbit anti-Human CD93 mAb
Remarks/Alias	C1QR1; C1qRP; CDw93; ECSM3; MXRA4; C1qR(P); dj737E23.1
Species	Human
Gene ID	22918
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-580 of human CD93 (NP_036204.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NPY3
Protein Weight	69kDa
Shipping	Ice bag