

Rabbit Anti-Human Cd94 Monoclonal Antibody For Fc - Q13241

Item Number: CM0011172



Introduction

Rabbit monoclonal antibody conjugated to , specific for human CD94 (KLRD1), validated for flow cytometry. Ideal for studying NK and T cell subsets, immune checkpoint regulation, and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD94 mAb
Remarks/Alias	KLRD1; CD94; natural killer cells antigen CD94
Species	Human
Gene ID (Human)	3824
Immunogen	Recombinant protein of human CD94.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13241
Protein Weight	17kDa/20kDa
Shipping	Ice bag