

Rabbit Anti-Human Cd94 Monoclonal Antibody For Flow Cytometry - Q13241

Item Number: CM0005894



Introduction

Rabbit monoclonal antibody against human CD94 (KLRD1), conjugated to , for flow cytometry. Recognizes the NK and T cell receptor CD94, involved in immune self-nonsel discrimination.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD94 mAb
Remarks/Alias	CD94
Species	Human
GeneID (Human)	3824
GeneID	3824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-179 of human CD94 (NP_002253.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13241
Protein Weight	21kDa
Shipping	Ice bag