

Rabbit Anti-Human Cd94 Monoclonal Antibody For Flow Cytometry - Q13241

Item Number: CM0004778



Introduction

-conjugated rabbit monoclonal antibody against human CD94, validated for flow cytometry. Ideal for NK and T cell immunophenotyping in immuno-oncology and infectious disease studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD94 mAb
Remarks/Alias	CD94
Species	Human
Gene ID	3824 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-179 of human CD94 (NP_002253.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	Q13241
Protein Weight	21 kDa
Shipping	Ice bag