

Rabbit Anti-Human Cd96/Tactile Monoclonal Antibody For Fc - P40200

Item Number: CM0009994



Introduction

-conjugated rabbit monoclonal antibody targeting human CD96 (TACTILE) for flow cytometry. Recombinant protein immunogen. Cross-reacts with human samples. Molecular weight: 65kDa. Ideal for studying T and NK cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD96/TACTILE mAb
Remarks/Alias	TACTILE
Species	Human
Gene ID	10225
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 25-519 of human CD96/TACTILE (NP_937839.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40200
Protein Weight	65kDa
Shipping	Ice bag