

Rabbit Anti-Human Cd96/Tactile Monoclonal Antibody For Flow Cytometry - P40200

Item Number: CM0008167



Introduction

-conjugated rabbit anti-human CD96/TACTILE monoclonal antibody, validated for flow cytometry. Recombinant immunogen (aa 25-519). Ideal for immune cell adhesion studies and NK/T-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD96/TACTILE mAb
Remarks/Alias	TACTILE
Species	Human
GeneID (human)	10225
GeneID	10225
Immunogen	Recombinant protein corresponding to a sequence within amino acids 25-519 of human CD96/TACTILE (NP_937839.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40200
Protein Weight	65kDa
Shipping	Ice bag