


Rabbit Anti-Human Cd97 Mab For Fc - P48960

Item Number: CM0006065



Introduction

Rabbit anti-Human CD97 monoclonal antibody for flow cytometry, reacts with human CD97. Useful for leukocyte adhesion and migration studies. Recognizes extracellular domain, 92kDa, SWISS P48960.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD97 mAb
Remarks/Alias	CD97; TM7LN1
Species	Human
GeneID (Human)	976
GeneID	976
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-398 of human CD97 (NP_510966.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P48960
Protein Weight	92kDa
Shipping	Ice bag