


Rabbit Anti-Human Cd97 Mab For Fc - P48960

Item Number: CM0005695



Introduction

-conjugated rabbit monoclonal antibody targeting human CD97 (ADGRE5) for flow cytometry. Involved in leukocyte adhesion and migration; broadly expressed on hematopoietic cells. Suitable for immunological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD97 mAb
Remarks/Alias	CD97; TM7LN1
Species	Human
GeneID (Human)	976
GeneID	976
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P48960
Protein Weight	81kDa/86kDa/92kDa
Shipping	Ice bag