

Rabbit Anti-Human Cd97 Monoclonal Antibody () For Flow Cytometry - P48960

Item Number: CM0012783



Introduction

Highly specific rabbit monoclonal antibody targeting human CD97 (ADGRE5). Conjugated to for direct immunofluorescence, ideal for flow cytometry analysis of hematopoietic cells and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD97 mAb
Remarks/Alias	ADGRE5; CD97; TM7LN1; CD97 antigen
Species	Human
Gene ID (Human)	976
Gene ID	976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-398 of human CD97 (NP_510966.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P48960
Protein Weight	81kDa/86kDa/91kDa
Shipping	Ice bag