

Rabbit Anti-Human Cd97 Monoclonal Antibody For Flow Cytometry - P48960

Item Number: CM0011452



Introduction

-conjugated rabbit monoclonal antibody targeting human CD97 (ADGRE5), a leukocyte adhesion GPCR. Validated for flow cytometry to study leukocyte migration, inflammation, and tumor biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD97 mAb
Remarks/Alias	CD97; TM7LN1
Species	Human
Gene ID (Human)	976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-398 of human CD97 (NP_510966.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P48960
Protein Weight	81kDa/86kDa/92kDa
Shipping	Ice bag