

Rabbit Anti-Human Cd98 Monoclonal Antibody For Fc - P08195

Item Number: CM0010357



Introduction

-conjugated rabbit monoclonal antibody targeting human CD98 (SLC3A2) for flow cytometry. Ideal for studying amino acid transport, integrin signaling, and viral entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD98 mAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Species	Human
Gene ID (Human)	6520
Gene ID	6520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-630 of human CD98 (NP_002385.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08195
Protein Weight	57kDa/61kDa/67kDa/68kDa/71kDa
Shipping	Ice bag