


Rabbit Anti-Human Cd99 Monoclonal Antibody For If/Icc, Fc - P14209

Item Number: CM0004963



Introduction

-conjugated rabbit monoclonal antibody against human CD99, validated for immunofluorescence (IF/ICC) and flow cytometry (FC). CD99 (MIC2) functions in T-cell adhesion and leukocyte transmigration. Suitable for human sample analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD99 mAb
Remarks/Alias	MIC2; HBA71; MIC2X; MIC2Y; MSK5X
Species	Human
Gene ID (Human)	4267
Gene ID	4267
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 24-123 of human CD99 (NP_002405.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P14209
Protein Weight	19kDa
Shipping	Ice bag