

Rabbit Anti-Human Cda Monoclonal Antibody For Fc (Intra) - P32320

Item Number: CM0010884



Introduction

-conjugated rabbit monoclonal antibody against human Cytidine deaminase (CDA/CDD) for intracellular flow cytometry. Reacts with human CDA, immunogen spans amino acids 1-146. Ideal for investigating nucleotide metabolism and granulocyte biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CDA mAb
Remarks/Alias	CDD
Species	Human
Gene ID (Human)	978
Gene ID	978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human CDA (NP_001776.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P32320
Protein Weight	16kDa
Shipping	Ice bag