

Rabbit Anti-Human Cda Monoclonal Antibody For Intracellular Flow Cytometry - P32320

Item Number: CM0010067



Introduction

-conjugated rabbit monoclonal antibody targeting human cytidine deaminase (CDA/CDD) for intracellular flow cytometry. Validated in human samples, it enables detection of CDA, a 16 kDa zinc-dependent enzyme highly expressed in granulocytes.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CDA mAb
Remarks/Alias	CDD
Species	Human
GeneID (Human)	978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human CDA (NP_001776.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P32320
Protein Weight	16kDa
Shipping	Ice bag