

Rabbit Anti-Human Cdcp1 Monoclonal Antibody For Fc - Q9H5V8

Item Number: CM0009209

Introduction

-conjugated rabbit monoclonal antibody targeting human CDCP1 (CD318, TRASK) for flow cytometry. Involved in cell adhesion, metastasis, and leukemia stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CDCP1 mAb
Remarks/Alias	CD318; TRASK; SIMA135
Species	Human
Gene ID (Human)	64866
Gene ID	64866
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9H5V8
Protein Weight	37kDa/72kDa/92kDa
Shipping	Ice bag