

Rabbit Anti-Human Ceacam5/Cd66E Mab For Fc - P06731

Item Number: CM0004601

Introduction

-conjugated rabbit monoclonal antibody targeting human CEACAM5/CD66e (CEA). Validated for flow cytometry (FC). Reacts with human. Immunogen: recombinant protein (aa 35-685). Ideal for cancer and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CEACAM5/CD66e mAb
Remarks/Alias	CEA; CD66e
Species	Human
GeneID (human)	1048
GeneID	1048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-685 of human CEACAM5/CD66e (NP_004354.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06731
Protein Weight	77kDa
Shipping	Ice bag