

Rabbit Anti-Human Ceacam5/Cd66E Monoclonal Antibody For Fc - P06731

Item Number: CM0004815



Introduction

Rabbit anti-Human CEACAM5/CD66e Monoclonal Antibody for flow cytometry. Recognizes human CEACAM5 (CEA/CD66e), a 77kDa cell surface glycoprotein, tumor marker, and cell adhesion molecule. Ideal for oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CEACAM5/CD66e mAb
Remarks/Alias	CEA; CD66e
Species	Human
GeneID (Human)	1048
GeneID	1048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-685 of human CEACAM5/CD66e (NP_004354.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06731
Protein Weight	77kDa
Shipping	Ice bag