

Rabbit Anti-Human Ceacam6 Monoclonal Antibody For Fc - P40199

Item Number: CM0009833



Introduction

-conjugated rabbit monoclonal antibody targeting human CEACAM6 (CD66c) for flow cytometry applications. CEACAM6 plays roles in cell adhesion and tumor progression; this antibody is suitable for oncology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CEACAM6 mAb
Remarks/Alias	NCA; CEAL; CD66c
Species	Human
Gene ID (Human)	4680
Gene ID	4680
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-320 of human CEACAM6 (NP_002474.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40199
Protein Weight	37kDa
Shipping	Ice bag