

Rabbit Anti-Human Ceacam6 Monoclonal Antibody For Flow Cytometry - P40199

Item Number: CM0004281



Introduction

-conjugated rabbit monoclonal antibody targeting human CEACAM6 (CD66c), a 37 kDa GPI-anchored cell adhesion molecule expressed in neutrophils and tumor cells. Suitable for flow cytometry (FC) applications.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CEACAM6 mAb
Remarks/Alias	NCA; CEAL; CD66c
Species	Human
Gene ID (Human)	4680
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-320 of human CEACAM6 (NP_002474.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40199
Protein Weight	37kDa
Shipping	Ice bag