

## Rabbit Anti-Human Claudin18.2 Monoclonal Antibody For Fc - P56856-2

Item Number: CM0008630

### Introduction

-conjugated rabbit monoclonal antibody targeting Claudin18.2 (isoform A2) validated for flow cytometry on human samples. Molecular weight 27kDa. Ideal for gastric tight junction and barrier function research.

[Learn More](#)

| Parameter        | Value                                                                                                                                         |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human Claudin18.2 mAb                                                                                                             |
| Remarks/Alias    | SFTA5; SFTPJ                                                                                                                                  |
| Species          | Human                                                                                                                                         |
| Gene ID (Human)  | 51208                                                                                                                                         |
| Gene ID          | 51208                                                                                                                                         |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human Claudin18.2(NP_001002026.1). |
| Source           | Rabbit                                                                                                                                        |
| Category         | Monoclonal Antibodies                                                                                                                         |
| Application      | FC                                                                                                                                            |
| Cross Reactivity | Human                                                                                                                                         |
| SWISS            | P56856-2                                                                                                                                      |
| Protein Weight   | 27kDa                                                                                                                                         |
| Shipping         | Ice bag                                                                                                                                       |