

Rabbit Anti-Human Claudin18.2 Monoclonal Antibody For Flow Cytometry - P56856-2

Item Number: CM0011314

Introduction

-conjugated rabbit monoclonal antibody targeting human Claudin18.2 (P56856-2). For flow cytometry. Immunogen: recombinant protein (aa 1-261). Reacts with human. Ideal for gastric barrier studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Claudin18.2 mAb
Remarks/Alias	SFTA5; SFTPJ
Species	Human
GenelD (human)	51208
GenelD	51208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human Claudin18.2 (NP_001002026.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P56856-2
Protein Weight	27kDa
Shipping	Ice bag