

Rabbit Anti-Human Clec12A/Cd371 Monoclonal Antibody For Fc - Q5Qgz9

Item Number: CM0008592



Introduction

-conjugated rabbit monoclonal antibody targeting human CLEC12A/CD371 (CLL-1), validated for flow cytometry. Myeloid inhibitory receptor involved in NETosis, inflammation, and pattern recognition. Ideal for immune profiling and functional studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CLEC12A/CD371 mAb
Remarks/Alias	CLEC12A; CD371; CLL-1; CLL1; DCAL-2; MICL; C-type lectin domain family 12 member A
Species	Human
GeneID (Human)	160364
GeneID	160364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-265 of human CLEC12A/CD371(NP_612210.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q5QGZ9
Protein Weight	8kDa/24kDa/27kDa/30kDa/32kDa
Shipping	Ice bag