

Rabbit Anti-Human Cpt1A Monoclonal Antibody For Flow Cytometry (Intracellular) - P50416

Item Number: CM0005963

Introduction

-conjugated rabbit monoclonal antibody targeting human CPT1A for intracellular flow cytometry. Detects 88 kDa liver isoform of carnitine palmitoyltransferase 1A, essential for fatty acid oxidation and mitochondrial metabolism.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CPT1A mAb
Remarks/Alias	CPT1; CPT1-L; L-CPT1
Species	Human
Gene ID (Human)	1374
Gene ID	1374
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 520-720 of human CPT1A (NP_001867.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P50416
Protein Weight	88kDa
Shipping	Ice bag