

Rabbit Anti-Human Cr2/Cd21 Monoclonal Antibody For Fc - P20023

Item Number: CM0008208



Introduction

-conjugated rabbit anti-human CR2/CD21 monoclonal antibody for flow cytometry. Targets complement receptor type 2 (SWISS P20023), a B-cell co-receptor and EBV receptor. Suitable for immune cell profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CR2/CD21 mAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9
Species	Human
GeneID (Human)	1380
GeneID	1380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-971 of human CD21 (NP_001868.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20023
Protein Weight	110kDa/113kDa/118kDa/119kDa
Shipping	Ice bag