

## Rabbit Anti-Human Cr2/Cd21 Monoclonal Antibody For Flow Cytometry - P20023

Item Number: CM0011984



### Introduction

-conjugated rabbit monoclonal antibody targeting human CR2/CD21, a B-cell complement receptor and EBV entry receptor. Validated for flow cytometry, it enables direct immunofluorescent detection of CD21-expressing cells in research applications.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CR2/CD21 mAb                                                                                   |
| Remarks/Alias    | CR; C3DR; CD21; CVID7; SLEB9                                                                                     |
| Species          | Human                                                                                                            |
| GeneID           | 1380                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 21-971 of human CD21 (NP_001868.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | FC                                                                                                               |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | P20023                                                                                                           |
| Protein Weight   | 110kDa/113kDa/118kDa/119kDa                                                                                      |
| Shipping         | Ice bag                                                                                                          |