

Rabbit Anti-Human Cytokeratin 18 Monoclonal Antibody For Intracellular Flow Cytometry (Fc) - P05783

Item Number: CM0011094

Introduction

-conjugated rabbit monoclonal antibody targeting human Cytokeratin 18 (CK-18, K18). Validated for intracellular flow cytometry. Suitable for studying epithelial cell biology, intermediate filaments, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 18 mAb
Remarks/Alias	K18; CK-18; CYK18
Species	Human
GeneID (Human)	3875
GeneID	3875
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 241-430 of human Cytokeratin 18 (NP_000215.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05783
Protein Weight	48kDa
Shipping	Ice bag