


Rabbit Anti-Human Cytokeratin 19 (Krt19) Monoclonal Antibody For Flow Cytometry - P08727

Item Number: CM0009002



Introduction

Rabbit anti-Human Cytokeratin 19 (KRT19) monoclonal antibody designed for intracellular flow cytometry applications. Suitable for epithelial cell biology and cancer stem cell research. - P08727

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 19 (KRT19) mAb
Remarks/Alias	K19; CK19; K1CS
Species	Human
GeneID (Human)	3880
GeneID	3880
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 301-400 of human Cytokeratin 19 (NP_002267.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P08727
Protein Weight	44kDa
Shipping	Ice bag