


Rabbit Anti-Human Cytokeratin 19 (Krt19) Monoclonal Antibody For Flow Cytometry - P08727

Item Number: CM0005921



Introduction

-conjugated rabbit monoclonal antibody targeting human Cytokeratin 19 (KRT19). Validated for intracellular flow cytometry. Suitable for epithelial biology, cancer research, and stem cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 19 (KRT19) mAb
Remarks/Alias	K19; CK19; K1CS
Species	Human
GeneID	3880
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 301-400 of human Cytokeratin 19 (KRT19) (NP_002267.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P08727
Protein Weight	44kDa
Shipping	Ice bag