


Rabbit Anti-Human Cytokeratin 8 Monoclonal Antibody For Flow Cytometry - P05787

Item Number: CM0010861



Introduction

-conjugated rabbit monoclonal antibody against human Cytokeratin 8 (K8, CK-8), a type II keratin intermediate filament protein, validated for intracellular flow cytometry (FC).

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 8 mAb
Remarks/Alias	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2
Species	Human
GeneID	3856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-483 of human Cytokeratin 8 (NP_002264.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05787
Protein Weight	53kDa/56kDa
Shipping	Ice bag