


Rabbit Anti-Human Cytokeratin 8 Monoclonal Antibody () For Intracellular Flow Cytometry - P05787

Item Number: CM0010459



Introduction

Rabbit monoclonal antibody against human Cytokeratin 8, conjugated to dye. Validated for intracellular flow cytometry (FC); targets KRT8, a key intermediate filament protein in epithelial cells. Ideal for cytoskeleton and epithelial biology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 8 mAb
Remarks/Alias	KRT8; CARD2; CK-8; CK8; CYK8; K2C8; K8; KO; keratin 8
Species	Human
GeneID (Human)	3856
GeneID	3856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-483 of human Cytokeratin 8 (NP_002264.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05787
Protein Weight	53kDa/56kDa
Shipping	Ice bag