


Rabbit Anti-Human Cytokeratin 8 Monoclonal Antibody For Intracellular Flow Cytometry - P05787

Item Number: CM0008285



Introduction

-conjugated rabbit anti-human cytokeratin 8 (CK8) monoclonal antibody for intracellular flow cytometry. Detects human KRT8, a simple epithelial marker, ideal for cancer and epithelial cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Cytokeratin 8 mAb
Remarks/Alias	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2
Species	Human
GeneID (Human)	3856
GeneID	3856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-483 of human Cytokeratin 8 (NP_002264.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05787
Protein Weight	53kDa/56kDa
Shipping	Ice bag