


Rabbit Anti-Human Dc-Signr/Cd299 Monoclonal Antibody For Flow Cytometry - Q9H2X3

Item Number: CM0012676



Introduction

Fluorescent -conjugated rabbit monoclonal antibody targeting human DC-SIGNR (CD299). Specific for flow cytometry (FC) applications, reactive with human. Used in viral attachment and immune surveillance research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human DC-SIGNR/CD299 mAb
Remarks/Alias	CD299; LSIGN; CD209L; L-SIGN; DCSIGNR; HP10347; DC-SIGN2; DC-SIGNR
Species	Human
GeneID (Human)	10332
GeneID	10332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 71-399 of human DC-SIGNR/CD299 (NP_055072.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9H2X3
Protein Weight	24kDa-45kDa
Shipping	Ice bag