

Rabbit Anti-Human Desmoglein 3/Pva Monoclonal Antibody For Flow Cytometry - P32926

Item Number: CM0010501



Introduction

-conjugated rabbit monoclonal antibody against human Desmoglein 3/PVA (DSG3), validated for flow cytometry. Targets the 130 kDa pemphigus vulgaris antigen, a desmosomal cadherin involved in cell adhesion. Ideal for skin and mucosal research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Desmoglein 3/PVA mAb
Remarks/Alias	PVA; ABOLM; CDHF6
Species	Human
Gene ID (Human)	1830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-615 of human Desmoglein 3/PVA (NP_001935.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32926
Protein Weight	108kDa
Shipping	Ice bag