


Rabbit Anti-Human Desmoglein 3/Pva Monoclonal Antibody For Flow Cytometry - P32926

Item Number: CM0004233



Introduction

High-performance -conjugated rabbit monoclonal antibody against human Desmoglein-3 (PVA/CDHF6). Validated for flow cytometry; ideal for investigating desmosomal adhesion, keratinocyte biology, and pemphigus vulgaris.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Desmoglein 3/PVA mAb
Remarks/Alias	PVA; ABOLM; CDHF6
Species	Human
Gene ID (Human)	1830
Gene ID	1830
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-615 of human Desmoglein 3/PVA (NP_001935.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32926
Protein Weight	108kDa
Shipping	Ice bag