

Rabbit Anti-Human Dnajb1/Hsp40 Monoclonal Antibody For Flow Cytometry - P25685

Item Number: CM0005936



Introduction

-conjugated rabbit monoclonal antibody targeting human DNAJB1/HSP40, validated for intracellular flow cytometry. Reacts with human samples; detects isoforms at 27 and 38 kDa. For research use in stress response and chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human DNAJB1/HSP40 mAb
Remarks/Alias	Hdj1; Sis1; HSPF1; Hsp40; RSPH16B
Species	Human
GeneID (Human)	3337
GeneID	3337
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 241-340 of human DNAJB1/HSP40 (NP_006136.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P25685
Protein Weight	27kDa/38kDa
Shipping	Ice bag