

Rabbit Anti-Human E-Cadherin/Cd324 Mab For Fc - P12830

Item Number: CM0012553



Introduction

ABflo 488-conjugated rabbit monoclonal antibody targeting human E-Cadherin/CD324, validated for flow cytometry (FC). Ideal for studying epithelial cell adhesion and cadherin-1 function. Cross-reacts with human. Predicted molecular weight 90/97 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human E-Cadherin/CD324 mAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324
Species	Human
GeneID (human)	999
GeneID	999
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 155-707 of human E-Cadherin/CD324 (NP_004351.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12830
Protein Weight	90kDa/97kDa
Shipping	Ice bag