

Rabbit Anti-Human E-Cadherin/Cd324 Monoclonal Antibody For Flow Cytometry (Fc) - P12830

Item Number: CM0010384



Introduction

-conjugated rabbit monoclonal antibody targeting human E-Cadherin/CD324 (SWISS: P12830). Validated for flow cytometry (FC), supporting research into cell adhesion, epithelial-mesenchymal transition, and tumor invasion.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human E-Cadherin/CD324 mAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324
Species	Human
Gene ID (Human)	999
Gene ID	999
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 155-707 of human E-Cadherin/CD324 mAb(NP_004351.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P12830
Protein Weight	97kDa
Shipping	Ice bag